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User Manual / Manual de instrucciones / Bedienungsanleitung /  
Manuale d'uso / Manual de instruções / Gebruiksaanwijzing

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---

## Thank you for choosing a Livoo product!

EN

Your blender, your barbecue, and your headphones are not just objects, they are extensions of yourself, allowing you to get closer to your friends and family, to bring them together, to laugh and to share experiences. Livoo's aim is to help you live your life to the fullest, whatever wishes you have throughout the year. Winter raclette, summer barbecue. Or the other way around. It's up to you. Take advantage of this freedom for feel good moments every day.

That's why we're constantly innovating and creating clever products full of zest, but most of all, accessible. Find our entire collection on our website [www.livoo.en](http://www.livoo.en).

## Livoo, your daily companion till social networks

@Livoo

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@Livoo

## Safety Instructions

Please, take time to read well the instruction sheet before the first use and keep them for a future use.

A misuse of the appliance can damage it or hurt the user. Assure that you use the appliance for the use for which it is designed. We refuse any responsibilities as for the damages by an incorrect use or a bad manipulation.

1. This appliance is intended for domestic use only. Any other use (medical, professional, commercial, etc.) is excluded.
2. Children should be supervised to ensure that they do not play with the appliance.
3. Keep this appliance indoors and in a dry place.
4. Do not immerse the unit in water or any other liquid.
5. Do not store any objects on the scale.
6. To avoid risk of injury, do not stand on the edges of the scale.
7. To prevent slipping, make sure that your feet and the scale are dry and free of dust before getting on the scale.
8. Do not disassemble the scale.
9. Avoid strong vibrations or shocks, such as dropping the scale on the floor.
10. Always place your scale on a stable, level surface. Avoid carpets and rugs.
11. Do not use the scale on a slippery surface or on a wet floor.
12. Keep the scale away from heat sources (radiators, etc.) to avoid deformation of plastic parts.
13. Do not store the scale in an upright position.
14. Use only the type of battery specified in this manual.
15. Replace the batteries when they are worn out or leak.
16. Remove the batteries if the appliance is not used for a long period of time.

# Technical Specifications

Batteries : 2 X AA included

In off mode, the device's consumption is 0W.

## Warning

For accurate weighing, we recommend weighing yourself every day at the same time, preferably without clothes and barefoot. Measurements of water or fat proportions, muscle or bone mass can only be done if you are barefoot.

This device is designed to determine body weight and analyze the proportion of water and fat in the body as well as muscle and bone mass. It can also calculate the body mass index (BMI).

Except for determining body weight, the other functions of the device cannot be used by individuals under 10 years old or over 70 years old.

The function for determining BMI is only suitable for adults aged 18 years and older. It will not work for anyone under 18 years old.

Please note that the results obtained from the scale are indicative. Different results may occur in individuals undergoing medical treatment or highly trained athletes. You should not base any treatment or diet on the results obtained from the scale. For that, you should consult your doctor.

The use of functions other than weighing is not recommended for the following people:

- People with a pacemaker or other medical implants.
- People showing symptoms of fever, edema, or osteoporosis.
- People undergoing dialysis.
- People being treated for cardiovascular diseases.

## Use of the device

### **Device Description (Figure 2)**

1	Handle
2	Thumb electrodes
3	Data cable
4	LED screen
5	Foot electrodes

Open the battery compartment cover on the back of the scale and insert the batteries while respecting the polarity. Replace the battery when "LO" appears on the screen, indicating that the battery is depleted. "----" appears on the screen to indicate an overload.

### **Weighing (without iOS/Android device) :**

Touch the scale with your foot to activate the screen. Wait until '0.0' (kg) appears. Note: For accurate results, please stand in the middle of the scale without moving. When the weight appears on the screen and flashes, you can step off the scale. The scale will turn off by itself after a few seconds.

### **Settings and storing personal data (Figure 1)**

To be able to use all the functions of your scale, you will need to download the "MovingLife" app from the App Store or Google Play Store. You can also scan the QR code to download the app (Figure 1).

### **Weighing with the app (Figures 3)**

Place the scale on a flat, hard surface. Step gently onto the scale barefoot (6) and place your feet on the 4 electrodes located on the left and right of the scale, then lift the scale's handle. Hold the handle with both hands (7), adjust the height of your hands (8), and straighten your arms (it is preferable that your arms do not touch your torso).

### Example of incorrect measurement posture (Figure 4)

A	Thumbs in contact
B	Both arms not stretched or extended
C	Measuring with socks
D	Squatting to measure
E	Ten fingers far from the electrodes
F	Knees bent or not straight
G	Both feet not placed correctly on the 4 electrodes

## Cleaning and maintenance

Never immerse the device in water. Do not use abrasive agents to clean it as this could damage the surface. Clean the surface of the scale with a soft damp cloth.

If you do not use the device for an extended period, remove the batteries and store it in a dry place, away from sunlight and out of reach of children.

The Bluetooth trademark belongs to Bluetooth SIG.

We

DELTA

BP61071

67452 Mundolsheim – France

#### **Declare under our own responsibility that the article :**

- **Description : Connected impedancemeter scale**
- **Reference : DOM531W**

to which this declaration applies, is in conformity with the applicable EU harmonization legislation :

DIRECTIVES	NORMES or SPECIFICATIONS
(RED) 2014/53/EU	ETSI EN 301 489-1 V2.2.3 (2019-11) ETSI EN 301 489-17 V3.2.4 (2020-09) ETSI EN 300 328 V2.2.2 (2019-07) EN IEC 62368-1 2020+A11 2020 EN 62479 2010 EN 50663 2021 EN 55014-1 2021 EN 55014-2 2021
ROHS 2011/65/EU	The object of the declaration described above is in conformity with Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment. EU 2011/65/EU et EU 2015/863

Mundolsheim, the 21/01/2025

David Bousquet, CEO

Signature:

### **Correct disposal of the unit (Electrical and Electronic Equipment).**

(Applicable in the European Union and the rest of Europe that have adopted separate collection systems countries). The European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal household waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain and reduce the impact on human health and the environment. The symbol of the crossed out dustbin is marked on all products to remind the obligations of separate collection. Consumers should contact their local authority or retailer for information concerning the correct disposal of their product.

- Non contractual illustration
- Product information is subject to change without prior notice.

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