
91916149/II-22
DE/EN/FR/ES/PT/IT/NL

Type 3750

Deutsch

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English

Product contains batteries and/or recyclable electric waste. For environment protection do not dispose in household waste, but for recycling take to electric waste collection points provided in your country.

Usage by children under age 3 is not recommended. Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge, if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and maintenance

shall not be performed by children. Children shall not play with the appliance. Keep away from children. Do not short-circuit battery terminals. Batteries may leak if empty or not used for a long time. In order to protect you and the appliance, please remove batteries in a timely manner and avoid skin contact whilst handling leaking batteries. Do not mix used and new batteries, different brands or types. Do not recharge non-rechargeable batteries. Exhausted batteries are to be removed from the appliance and safely disposed of.

PRIOR TO USE, ENSURE THE BRUSH HEAD FITS WELL. REPLACE BRUSH HEAD EVERY 3 MONTHS OR SOONER IF BRUSH HEAD BECOMES WORN OR NO LONGER FITS PROPERLY OR IF THE BRUSH IS DROPPED.

Do not use as a toy or without adult supervision. Small parts may pose

choking hazard. Do not chew or bite on brush. For cleaning, switch off the handle and remove the brush head. Rinse both parts separately with water and wipe them dry before reassembling.

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